

Putting Insets into Plots

VNMRj 2.1A walkup mode

1. Load the file and process the spectrum.
2. Plot the whole spectrum or large expansion by using the **Process->Plot** menu such as [**Plot Spectrum**] [**Plot Scale**]. Do not hit the [**Plot Page**] button yet.
3. Go to the **Process->Line Lists** menu. Put the red cursors around the peak(s) for the inset. Select [**Inset Spectrum**].
4. The peak for the inset will be displayed above the spectrum. The vertical scale of the inset can be changed with the middle mouse button at this point. To change the position of the inset and to expand it, select the shift position button (blue box) from the icon menu on the left. After clicking the shift position button, the left mouse button will move the inset right or left and the right mouse button will enlarge the inset.



5. After adjusting the inset, select [**Plot inset**]. Then go back to the **Process->Plot** page and select [**Plot Page**] to finish the plot.

